



[Altogen Biosystems](#)

848 Rainbow Blvd Suite 823

Las Vegas, NV 89107 USA

Tel: (702) 349-6103

Fax: (702) 989-0841

Web: www.altogen.com

E-mail: orders@altogen.com

Lungs *In Vivo* Transfection Reagent

Cat. No. 5080	Size: 0.5 ml	Cat. No. 5082	Size: 8.0 ml
Cat. No. 5081	Size: 1.5 ml	Cat. No. 5083	Size: 25 ml

Contents and Shipping:

Lungs *in vivo* transfection reagent is supplied in liquid form at a concentration of 1.5 mg/ml. Shipped at ambient temperature. Transfection Enhancer reagent is supplied in liquid form at a concentration of 1 mg/ml.

Description:

Lungs transfection reagent is a nanoparticle liposome-based formulation optimized for *in vivo* administration with targeted biomolecules delivery into lungs and lung tumors. Transfection reagent is a proprietary, animal-origin-free formulation with high transfection efficiency of small RNA (siRNA, shRNA, microRNA) and plasmid DNA.

Storage:

Store reagent at 4°C upon receipt. If stored properly, reagent is stable for 6 months.

Product Qualification:

Reagents are tested for functional activity, absence of nuclease contamination and microbial contamination.

Reagent Properties:

Transfection complexes are stable in serum for at least 8 hours.

Efficient delivery to the lung tissue and lung tumors via intranasal and inhalation administration.

Minimal toxicity.

Recommended Modes of Administration:

Intranasal administration.

Inhalation administration.

Intended Use:

For research use only. Not intended for any diagnostic or therapeutic use.

MSDS:

MSDS documents are available online at www.altogen.com

General Guidelines:

Handle animals and conduct experiments according to national regulations and approvals by the local IACUC and ethics committee. Please note that mouse and rat administrations may require skills, experience, and associated equipment (e.g. restraining device). All persons handling animals should be properly trained.

Guidelines for Small RNA (siRNA, shRNA, microRNA) Administration:

Prepare siRNA stock solution in DNase- and RNase-free water.

Recommended dose of siRNA (per injection):

- 1-2 µg siRNA with chemical modification (resistant to nuclease degradation *in vivo*), or
- 25-30 µg standard, non-modified siRNA

Recommended Protocol (for 1 injection):

1. Dilute 30 µg of plasmid DNA and/or 1-2 µg of siRNA (chemically modified) in 35 µl RNase/DNase-free water. Vortex gently.
2. Add 35 µl of diluted DNA (and/or RNA) to the sterile tube containing 35 µl Transfection Reagent.
3. Incubate for 10-15 minutes at room temperature.
4. Add 10 µl of Transfection Enhancer Reagent. Vortex gently.
5. Incubate for 2-3 minutes at room temperature.
6. Perform direct intranasal administration (40 µl per each mouse nostril), or add required amount of sterile solution of 5% glucose (w/v) to perform inhalation administration (please see individual inhaler protocol for recommended settings typically available in the inhaler manual).
7. Please note that delivery efficiency can be significantly increased by performing secondary administration (at least 8-12 hours after first administration).
8. Assay for target gene expression. Maximum mRNA expression effect is typically observed 8-24 hours after administration, while maximum effect on protein level is achieved 18-36 hours post-administration.

Dosing:

Animal body weight (g)	Transfection Reagent (µl) / per single administration
10 – 20 g (mouse)	35 µl
20 – 50 g (small rodents)	70 µl
50 – 100 g (small rodents)	105 µl
150 – 200 g (small rodents)	140 µl
200 – 250 g (small rodents)	175 µl
250 – 300 g (rat)	210 µl
300 – 350 g (rat)	245 µl
350 – 400 g (rat)	280 µl

Limited Use Label License:

The purchase of this product conveys to the purchaser the limited right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. This product is for research purposes only and is not for use in commercial applications of any kind. For information on obtaining additional rights, please contact Altogen Biosystems at orders@altogen.com

Limited Product Warranty:

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Altogen Biosystems. Altogen Biosystems shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product. This product is developed, manufactured and sold for research purposes only. Reagent is not suitable for administration to humans.